



Recirculating heat and many other exclusive features



HEPA filters



Certified class 100 High Efficiency. Particulate air (HEPA) filters remove 99.97% of all particles larger than 0.3 microns that transport bacteria and viruses. Individually tested to assure entire HEPA filter system is leak-free.

Removable trays



Model 44 has nine large 18" x 18" (46 x 46 cm) parts trays standard. Model 43 has three 18" x 18" (46 x 46 cm) trays and six 18" x 9" (46 x 23 cm) half trays standard.

Bag Drier



Optional Bag Drier solves the problem of drying both the inside and outside of breathing bags after pasteurization. Slides into any Sterile-Drier model manufactured since 1984 and dries up to 16 bags at once.

Slide-out tube holder



Easy-loading tube-holders slide out at convenient height. Holds 35 tubes up to 5 ft. (1.5 m) long. Tube-holders can be replaced with more parts trays or Bag Drier.

Features and Specifications

Safety features

Over-temperature thermal cut-offs on heaters and blower motors. Individually fused, 3-wire hospital-grade electrical connectors. Each Sterile-Drier individually tested for 1350-volt dielectric insulation

Dimensions

Model 43: 26"W x 20.5"D x 86"H (66 x 52 x 218 cm)
Model 44: 46.5"W x 20.5"D x 86"H (118 x 52 x 218 cm)

Ordering Information

Olympic Sterile-Driers

Model 43: With 3 full-size parts drying trays, 6 halfsize parts drying trays, and tube-holder for 35 tubes. 120V~, 60 Hz, 900 watts
CAT. NO. 54343

Model 44: With 9 full-size parts drying trays and tubeholder for 35 tubes. 120V~, 60 Hz, 1700 watts
CAT. NO. 54344

Note: Sterile-Drier models are available with optional rear doors and additional trays and tube holders

Accessories

Bag Drier: For use with any Olympic Sterile-Drier manufactured since 1984
CAT. NO. 54355

natus

Natus Medical Incorporated
1501 Industrial Road
San Carlos, CA 94070 USA
1-800-303-0306
+1-650-802-0400
www.natus.com